

MACHACADORA DE MANDÍBULAS JC-50

Instrumento de laboratorio para la molienda seca de material de dureza media-alta como granito, feldspatos, chamota, gres... etc.
Producción de polvo de granulometría media-fina.
Fácil intercambio de mandíbulas.

Especificaciones técnicas:

Máxima dureza del material a moler: 6 – 7 escala MOHS
Tolva de carga 180x190mm con alimentador de 170mm
Rueda manual para la regulación de la machacadora
Producción: 10 – 50Kg/h
Tamaño máximo de muestra: 5 – 20mm
Alimentación: 400V – 50Hz trifásico

Código GT0582

Modelo JC-50

Dimensiones exteriores 800x850x680

Potencia 0,75KW

Peso 205Kg



JAW CRUSHER JC-50

Laboratory instrument for the dry grinding of medium– high hardness material such as granite, feldspar, chamotte, gres, etc.

Production of grit with medium-fine granulometry .Easy change hardened steel jaws

Technical specifications:

Max hardness grinding's material: 6 ÷ 7 MOHS scale.

Loading hopper 180x190 mm with feed of 170 mm.

Hand wheel for crusher regulation.

Production: 10 ÷ 50 kg/h.

Feeding max sample size: 5 ÷ 20 mm.

Supply: 400 V - 50 Hz three phases

Code GT0582

Model JC-50

Overall dimensions: 800x850x680

Power 0,75KW

Weight 205Kg