

Surface Tension

*– rapid, easy,
and cost-effective*

AquaPi



AquaPi

Ideal for laboratory or on-site measurements in many industries.
Extremely simple and easy to use.




- Fast, automated measurement
- Handy for basic research in surface chemistry
- Ideal for checking washing process results
- Useful for water quality and environmental monitoring
- Simple tensiometer for educational purposes
- Single-use sample cups reduce contamination risk

Specifications

Sample cup volume	3 ml
Measurement range	10 – 100 mN/m
Accuracy/sensitivity	better than 0.1 mN/m
Dimensions	30 × 8 × 15 cm (11.8 × 3.2 × 5.9 in)
Weight	2.2 kg (4,9 lb)

International distribution



Kibron Inc.
P. O. Box 141
FIN-00171 Helsinki FINLAND
Tel: +358 9 6811 920
Fax: +358 9 6811 9222
Email: info@kibron.com
www.kibron.com