



## MOLINO DE ASPAS MPT

Código / Code	Modelo / Model	Dimensiones exteriores / Overall dimensions	Dimensiones cámara de molienda / Crushing Chamber dimensions	RPM	Potencia / Power	Peso / Weight
310.0	MPT0	150x150x300	78x113 (150g)	30000	820W	5Kg
310.1	MPT1	180x160x340	85x139 (250g)	24000	1200W	8kG
310.2	MPT2	200x170x390	140x174 (350g)	24000	1800W	12Kg