

ramé-hart Model 190

Contact Angle (CA) Goniometer with DROPimage CA Software

ramé-hart Model 190 (p/n 190-F2)

Model 190 is an affordable and easy-to-use contact angle tool designed to offer the highest level of precision computer-aided contact angle analysis at the lowest possible acquisition cost. Model 190 accepts the Manual Tilting Base, Environmental Fixture, Wafer Support up to 4" in diameter, and other options. This turn-key system includes DROPimage CA software.

What's in the box: Goniometer with F2 Series Camera, Microsyringe Fixture and Shade, LED Backlight Illuminator, power supply, (1) Microsyringe Assembly and (1) 22g Straight Needle for manual dispensing, PC with FireWire Interface and LCD, DROPimage CA Software, User Manual, Calibration Tool, and Storage Cover — everything needed to start taking contact angle measurements right out of the box.



Specifications

Stage Size 2 x 3 in (51 x 76mm)

Sample Size 3 in (76mm) deep x unlimited

Contact Angle Range 0 to 180°

Resolution 0.1°

Accuracy +/- 0.10°

Camera IEEE1394a FireWire, 60 fps

Backlighting LED Illuminator

Stage Precision 3-Axis Locking with

Leveling Specimen Stage

Dimensions 19 x 20 x 10 in (480x500x250mm)

Weight 14 lbs / 6.3 kg (excludes PC)

Power Supply 110 or 220 VAC

DROPimage Features Contact Angle Measurement

Calibration Tool

Spares Microsyringe Assy. 100-10-20

Needle 100-10-12-22

Options Environmental Fixture 100-07-60

Wafer Support (4") 100-21-4 Software Upgrade 400-22-30

Manual Tilting Base 100-25-M

Film Clamps 100-15

Note: The Automated Dispensing System (100-22) can be added to Model 190 but requires the DROPimage Standard or Advanced Software Upgrade.

ramé-hart instrument co · www.ramehart.com · carl@ramehart.com · 973-448-0305 · fax 0315

19 Route 10 East · Suite 11 · Succasunna · New Jersey · 07876 · USA

